



COLGATE-PALMOLIVE Company
MATERIAL SAFETY DATA SHEET

This industrial Material Safety Data Sheet is not intended for consumers and does not address consumer use of the product. For information regarding consumer applications of this product, refer to the product label.

**Colgate Visible White Chairside Procedure Kit
 Post -Treatment Gel**

Effective Date: June 14, 2007

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

<p>COLGATE-Oral Pharmaceuticals, INC. 300 Park Avenue New York NY 10022</p> <p>PRODUCT NAME: Colgate Visible White Chairside Procedure Kit /Post -Treatment Gel CAS NUMBER: Not applicable – product is a mixture GENERAL USE: A formulated post whitening treatment gel.</p>	<p>EMERGENCY TELEPHONE NUMBER: For emergency involving spill, leak, fire, exposure or accident, call CHEMTREC (800) 424-9300, day or night</p> <p>FOR MEDICAL EMERGENCIES INVOLVING THIS PRODUCT CALL (888) 489-3861</p>
--	--

<p>2. COMPOSITION/INFORMATION ON INGREDIENTS</p> <p>OSHA-REGULATED COMPONENTS (present at a concentration of > or = 1%):</p> <table border="1"> <thead> <tr> <th>Component</th> <th>CAS</th> <th>%</th> <th>PEL</th> <th>TLV</th> </tr> </thead> <tbody> <tr> <td>Sodium Hydroxide</td> <td>1310-73-2</td> <td>5-10</td> <td>2 mg/m³</td> <td>2 mg/m³</td> </tr> </tbody> </table>	Component	CAS	%	PEL	TLV	Sodium Hydroxide	1310-73-2	5-10	2 mg/m ³	2 mg/m ³	<p>The following components, present at a concentration of > or = 0.1%, are listed as carcinogens or potential carcinogens by either the National Toxicology Program (NTP), the International Agency for Research on Cancer (IARC) or OSHA:</p> <table border="1"> <thead> <tr> <th>Component</th> <th>CAS#</th> <th>%</th> <th>PEL</th> <th>TLV</th> </tr> </thead> <tbody> <tr> <td>Not Applicable</td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Component	CAS#	%	PEL	TLV	Not Applicable				
Component	CAS	%	PEL	TLV																	
Sodium Hydroxide	1310-73-2	5-10	2 mg/m ³	2 mg/m ³																	
Component	CAS#	%	PEL	TLV																	
Not Applicable																					

3. HAZARDS IDENTIFICATION

POTENTIAL HEALTH EFFECTS

EYE CONTACT: Causes severe eye irritation on direct contact.

SKIN CONTACT: May cause skin irritation on prolonged or repeated contact.

INGESTION: May be harmful if swallowed in large amounts.

INHALATION: No adverse effects expected.

4. FIRST AID MEASURES

EYE CONTACT: Flush eyes with large amounts of water for 15 minutes. Get medical attention if irritation persists.

SKIN CONTACT: Rinse area with plenty of water. Get medical attention if irritation persists.

INGESTION: Drink 1-2 glasses of water. Get medical attention.

INHALATION: Not Applicable.

5. FIRE FIGHTING MEASURES

EXTINGUISHING MEDIA: Water spray, all-purpose dry chemical, CO₂.

SPECIAL FIRE FIGHTING PROCEDURES: Self-contained breathing apparatus and protective clothing should be worn when fighting chemical fires.

6. ACCIDENTAL RELEASE MEASURES

Cover with inert, absorbent material and remove to disposal container. Spill area may be slippery. Flush with plenty of water.

7. HANDLING AND STORAGE

Store in a tightly closed container in a cool, dry, well-ventilated area.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

GENERAL CONTROLS: Avoid eye contact.

PROTECTIVE CLOTHING: The use of safety goggles and protective gloves is recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance and odor: Clear, colorless gel with no characteristic odor.

pH: ~6.5 – 7.5

Flash Point (Method): Not Applicable

10. STABILITY AND REACTIVITY

GENERAL: This product is stable. Hazardous polymerization will not occur.

INCOMPATIBLE MATERIALS: Avoid contact with strong caustics and most metals.

HAZARDOUS DECOMPOSITION: None known.

11. TOXICOLOGICAL INFORMATION

This product has not been tested as a whole. However, this formula was reviewed by expert toxicologists in the Product Safety Assurance Department of Colgate-Palmolive and is determined to be safe for its intended use. This review has taken into consideration available safety-related information including information on individual ingredients, similar ingredients, similar formulas and potential ingredient interactions. This review is a component of the hazard determination used to prepare the statements in Section 3 of the MSDS.

12. ECOLOGICAL INFORMATION

There is currently no ecological information available at this time.

13. DISPOSAL CONSIDERATIONS

Any disposal practice must be in compliance with local, state, and federal laws and regulations (contact local or state environmental agency for specific rules). Do not dump into sewers, any body of water or onto the ground.

14. TRANSPORTATION

DOT designation: UN3139, Oxidizing liquids, n.o.s. (contains Potassium nitrate), 5.1, PG III

15. REGULATORY INFORMATION

<p>RCRA (40 CFR 261, Subpart D): Not Applicable</p> <p>CLEAN WATER ACT: Contains sodium fluoride and sodium hydroxide which are Section 311 materials.</p> <p>CLEAN AIR ACT: Contains polyethylene glycol which is a Section 111 material.</p> <p>SARA: Sections 301-304 (Threshold planning quantity - TPQ) 40 CFR 355: Not Applicable</p> <p>Section 313 (Toxic chemical release reporting) 40 CFR 372: The following chemicals must be reported under SARA 313: Not Applicable</p> <p>CERCLA: Section 102 (Reportable Quantity - RQ) 40 CFR 302: The Reportable Quantity (RQ) for this product to the environment is 18,519 lbs. based on the presence of sodium hydroxide (~5.4%). Releases greater than or equal to 18,519 lbs. must be reported to the National Response Center (NRC) immediately: 800-424-8802.</p>	<p>New Jersey Right to Know Hazardous Substance List: This product contains the following components subject to reporting requirements: Sodium hydroxide, potassium nitrate, sodium fluoride</p> <p>Pennsylvania Hazardous Substance List: This product contains the following components subject to reporting requirements: Glycerin, sodium hydroxide, potassium nitrate, sodium fluoride</p> <p>Massachusetts Substance List: This product contains the following components subject to reporting requirements: Sodium hydroxide, potassium nitrate, sodium fluoride</p> <p>Canada: Workplace Hazardous materials Information System (WHMIS)-listed material. This product contains the following components subject to reporting requirements: Sodium hydroxide</p>
--	---

16. OTHER INFORMATION

Effective Date: June 14, 2007

MSDS Status: New

The information on this sheet is limited to the material identified and is believed by the Colgate-Palmolive Company to be correct based on its knowledge and information as of the date noted. Colgate makes no representation, guarantee or warranty, expressed or implied, as to the accuracy, reliability or completeness of the information and assumes no responsibility for injury, damage or loss resulting from the use of the material.